

## GUIDELINES FOR SEQUENCING SAMPLE SUBMISSION


**DNA Type - Genomic DNA; BAC, Plasmid, Cosmid DNA; Long-range PCR (>1.5kb)**

DNA Quantity for library preparation:

- At least 2ug for Tru-Seq Illumina DNA library preparation
- Quantitation should be performed using a **PicoGreen** Assay (if you do not have one, we can do it for you)
- ***NanoDrop UV and NanoDrop Fluorometry are NOT compatible with this protocol***
- We highly recommend you to prepare and provide us some extra amount of DNA to ensure timely process of your samples

DNA Quality for library preparation:

- Must be double-stranded
- Should not be a product of whole genome amplification
- Should not be degraded – verified by agarose gel
- ***The DNA sample will NOT be accepted with out the agarose gel image***
- Should contain NO particulate matter
- Should have an OD 260/280 ~ 1.8
- Sample should have a minimal concentration of 50ng/ul in TE, H<sub>2</sub>O, or Tris.
- Sample should NOT have RNA contamination (gDNA or BAC DNA) or primer/primer dimers (long-range PCR)



Example of gDNA run in 1% agarose gel

## RNA Type

### The sample RNA Requirements for Tru-Seq Illumina RNA library preparation:

- total amount of RNA 2ug
- quantitated by Ribogreen
- Sample should have a minimal concentration of 50ng/ul in H2O
- pure (OD 260/280 > 2.0)
- DNA free
- quality assessed on an RNA 6000 Pico Chip on the Agilent 2100 Bioanalyzer instrument. A typical RNA sample will produce a smear that ranges from 0.2 kb to 7 kb. This protocol is not designed for preparing small RNA molecules, for example snoRNA, microRNA, tRNA, etc.

For all other types of samples, please contact Alina Akhunova via email [akhunova@ksu.edu](mailto:akhunova@ksu.edu) or by phone (785) 532-1393